

TeSR™-E8™

Xeno-Free and Feeder-Free Culture Medium for hPSCs

TeSR™-E8™ is a highly defined, xeno-free and feeder-free culture medium for human pluripotent stem cells (hPSCs). It is based on the E8 formulation (Table 1) published by the laboratory of Dr. James Thomson (University of Wisconsin-Madison), the lead research group behind the design of mTeSR™1.¹⁻³

Like mTeSR™1, TeSR™-E8™ is made with the highest level of quality and care. Xeno-free TeSR™-E8™ contains only the essential components required for maintenance of embryonic stem (ES) and induced pluripotent stem (iPS) cells, thus providing a more simple medium for hPSC culture. It can be used with a surface coating of Corning® Matrigel® hESC-Qualified Matrix or with Vitronectin XF™.

Cells cultured using TeSR™-E8™ maintain typical pluripotent stem cell morphology (Figure 1), display high expansion rates (Figure 2), express high levels of undifferentiated cell markers (Figure 3), exhibit normal karyotype after long-term passaging (www.stemcell.com/TeSR-E8Data) and retain the ability to differentiate into cells of all three germ layers (Figure 4).




FIGURE 1. Morphology of Human ES and iPS Cells in TeSR™-E8™ Cultures

Undifferentiated human **(A)** ES (H9) and **(B)** iPS (WLS-1C) cells cultured on Corning® Matrigel® in TeSR™-E8™ retain the prominent nucleoli and high nuclear to cytoplasm ratio characteristic of this cell type. Densely packed cells and multilayering are prominent when cells are ready to passage.

TeSR™-E8™ is part of the most complete, defined system of reagents for hPSC research:

- **mTeSR™1, TeSR™2 and TeSR™-E8™** media for feeder-free culture
- **ReLeSR™** for easy passaging of hPSCs without manual selection or scraping
- **Vitronectin XF™** for a completely xeno-free culture system when used with TeSR™2 or TeSR™-E8™
- **STEMdiff™** suite of products for optimized hPSC differentiation to specific lineages
- **CryoStor®** or **FreSR™** media for cryopreservation




FIGURE 2. High Expansion Rates are Observed with TeSR™-E8™

Graph shows the average fold expansion per passage ± SEM obtained for human ES and iPS cells cultured in TeSR™-E8™ with Corning® Matrigel®, and passaged using ReLeSR™ or Gentle Cell Dissociation Reagent (GCDR) protocols over 10 passages. Expansion was determined by enumerating the cell aggregates obtained at harvest and dividing by the number of cell aggregates seeded.

Note: This data is representative of cultures passaged after 6-7 days in culture; lower expansion should be expected if using shorter culture times.

TeSR™-E8™

Feeder-Free Culture Medium for Human Pluripotent Stem Cells




FIGURE 3. Cells Cultured in TeSR™-E8™ Medium Express Markers of Undifferentiated Pluripotent Stem Cells

(A) Histogram analysis for H9 human ES and WLS-1C human iPS cells characterized using flow cytometry for markers of undifferentiated cells (SSEA-3 and OCT4) after 10 passages in TeSR™-E8™ (blue or red = sample, black = secondary antibody only). **(B)** H1 cells cultured in TeSR™-E8™ medium were characterized by OCT4 (OCT3) immunostaining after 18 passages.

Product Information

PRODUCT	CATALOG #
TeSR™-E8™	05940
ReLeSR™	05872
Vitronectin XF™	07180; 07190
Non-Tissue Culture Treated 6-Well Plates	27147
mFreSR™	05855
CryoStor® CS10	07930
CryoStor® CS5	07933
CryoStor® CS2	07932

References

- Ludwig TE et al. (2006) Nat Methods 3(8): 637–46.
- Ludwig TE et al. (2006) Nat Biotechnol 24(2): 185–7.
- Chen G et al. (2011) Nat Methods 8(5): 424–9.

Copyright © 2016 by STEMCELL Technologies Inc. All rights reserved including graphics and images. STEMCELL Technologies & Design, STEMCELL Shield Design, Scientists Helping Scientists, FreSR, ReLeSR and STEMdiff are trademarks of STEMCELL Technologies Canada Inc. TeSR, E8 and mTeSR are trademarks of WARF. Corning and Matrigel are trademarks of Corning Inc. Vitronectin XF is developed and manufactured by Primorigen Biosciences Inc., and Vitronectin XF™ is a trademark of Primorigen Biosciences, Inc. CryoStor is a registered trademark of BioLife Solutions. All other trademarks are the property of their respective holders. While STEMCELL has made all reasonable efforts to ensure that the information provided by STEMCELL and its suppliers is correct, it makes no warranties or representations as to the accuracy or completeness of such information.

STEMCELL TECHNOLOGIES INC.'S QUALITY MANAGEMENT SYSTEM IS CERTIFIED TO ISO 13485. PRODUCTS ARE FOR RESEARCH USE ONLY AND NOT INTENDED FOR HUMAN OR ANIMAL DIAGNOSTIC OR THERAPEUTIC USES UNLESS OTHERWISE STATED.




FIGURE 4. Directed Differentiation of Human iPS Cells Cultured in TeSR™-E8™

STiPS-F016 human iPS cells maintained in TeSR™-E8™ were differentiated into all three germ layers. Endoderm specification was achieved using the STEMdiff™ Definitive Endoderm Kit (TeSR™-E8™ Optimized) (Catalog #05115). Mesoderm specification was demonstrated using STEMdiff™ Mesoderm Induction Medium (Catalog #05220). Cells were analyzed by flow cytometry for expression of endoderm (CXCR4 and SOX17) or mesoderm (Brachyury) markers. Ectoderm specification was demonstrated using STEMdiff™ Neural Induction Medium (Catalog #05835). Central Nervous System (CNS)-enriched neural progenitor cell (NPC) cultures expressing PAX6 (green) are shown.

COMPONENT	mTeSR™1	E8*
DMEM/F12 (DF12), NaHCO ₃ , L-Ascorbic Acid, Selenium, Transferrin, Insulin, FGF2, TGF-β	•	•
Bovine Serum Albumin (BSA)	•	
Glutathione	•	
β-mercaptoethanol (BME)	•	
Pipecolic Acid	•	
GABA	•	
Lithium Chloride	•	
Defined Lipids	•	
Trace Elements	•	

TABLE 1. TeSR™-E8™ Product Specifications

*As published in Chen et al. 2011, the base design for TeSR™-E8™.

For a complete list of related products, including specialized cell culture and storage media, matrices, antibodies, cytokines and small molecules, visit www.stemcell.com/hPSCworkflow or contact us at techsupport@stemcell.com. Recommended antibodies include Anti-Human OCT4 (OCT3) Antibody, Clone 3A2A20 (Catalog #60093).