

FETAL BOVINE SERUM

**FOR USE IN METHYLCELLULOSE ASSAYS OF MURINE
MYELOID COLONY-FORMING CELLS**

Catalog # 06200	100 mL per bottle
Catalog # 06250	500 mL per bottle

PRODUCT DESCRIPTION:

This fetal bovine serum has been tested with all of our culture reagents used in the preparation of methylcellulose assays of murine myeloid colony-forming cells (CFU-E, BFU-E, CFU-G/M and CFU-GEMM).

Components combined with this product should be pre-tested in your methylcellulose culture system and selected to provide satisfactory colony growth.

All sera have been pre-filtered.

STABILITY/STORAGE:

Store at -20°C.

Stable for two months when stored at 4°C.

Contents guaranteed sterile if seal is not tampered with.

**THIS REAGENT IS FOR RESEARCH USE ONLY.
IT IS NOT TO BE ADMINISTERED TO HUMANS.**