


PRODUCT DESCRIPTION:

This product has been developed for use in applications where a culture medium of defined composition, containing only human-derived components, is required. This product contains pre-screened batches of human serum albumin (HSA) that have been selected to support the expansion of human hematopoietic progenitor cells in serum-free media formulations.

Product Information Sheet



**HUMAN
SERUM
ALBUMIN**

**10% HSA in
Iscove's MDM**

Version 1.0.0

Catalog #09350 100 mL

FORMULATION:

Contains 10% HSA in Iscove's MDM.

This product contains biological reagents, and as such cannot be completely characterized or quantified. Some variability is unavoidable.

This product is aseptically manufactured using tightly controlled processes and pre-screened components, and is sterility tested.

STABILITY AND STORAGE:

Product stable at -20°C until expiry date indicated on label. Stable for 2 months at 2-8°C.

It is recommended that the product be thawed, mixed well and aliquoted into smaller volumes. These aliquots should be stored at -20°C until use. *Avoid freezing and thawing repeatedly.*

DIRECTIONS FOR USE:

Use at a final HSA concentration of 1 - 2% in media containing desired supplements, for culture of human hematopoietic progenitors.

See Material Safety Data Sheet for more information.

**FOR RESEARCH USE ONLY.
NOT FOR THERAPEUTIC OR DIAGNOSTIC USE.**

StemCell Technologies

In North America

Tel: 1.604.877.0713
Fax: 1.604.877.0704
Toll Free Tel: 1.800.667.0322
Toll Free Fax: 1.800.567.2899
e-mail: info@stemcell.com
www.stemcell.com

In the United Kingdom

Tel: +44.(0).20.7691.3561
Fax: +33.(0).4.76.18.99.63
Toll Free within United Kingdom:
Tel: 0800.731.27.14
Fax: 0800.731.27.13
e-mail: info@stemcellgb.com

In Europe

Tel: +33.(0).4.76.04.75.30
Fax: +33.(0).4.76.18.99.63
e-mail: info@stemcellfrance.com

*Revised:
March 2008*