

## Release Notes for STEMvision™ Software Version 2016.04.21.00

These release notes describe the new functionality and/or fixes to the previous STEMvision™ Software Version 2015.12.18.01

SOFTWARE	NEW FUNCTIONALITY/FIX
STEMvision™ Acquisition	Version 2.0.1.0
STEMvision™ Analyzer (Human)	Version 2.0.2.2 <ul style="list-style-type: none"> <li>• Sending end of report data to HL7 server</li> </ul>
STEMvision™ Colony Marker	Version 2.0.3.0 <ul style="list-style-type: none"> <li>• Fix a bug that is not creating a summary image for manually marked images (colour and mono).</li> <li>• Fix a bug that is not using a correct image for manually marked images.</li> </ul>

Copyright © 2016 by STEMCELL Technologies Inc. All rights reserved including graphics and images. STEMCELL Technologies & Design, STEMCELL Shield Design, Scientists Helping Scientists, and STEMvision are trademarks of STEMCELL Technologies Inc. While STEMCELL has made all reasonable efforts to ensure that the information provided by STEMCELL and its suppliers is correct, it makes no warranties or representations as to the accuracy or completeness of such information.

STEMCELL TECHNOLOGIES INC.'S QUALITY MANAGEMENT SYSTEM IS CERTIFIED TO ISO 13485. PRODUCTS ARE FOR RESEARCH USE ONLY AND NOT INTENDED FOR HUMAN OR ANIMAL DIAGNOSTIC OR THERAPEUTIC USES UNLESS OTHERWISE STATED.